



# PDA100A2

## Si Switchable Gain Detector

### User Guide



## Table of Contents

<b>Chapter 1</b>	<b>Warning Symbol Definitions .....</b>	<b>1</b>
<b>Chapter 2</b>	<b>Description.....</b>	<b>2</b>
<b>Chapter 3</b>	<b>Setup.....</b>	<b>3</b>
<b>Chapter 4</b>	<b>Operation.....</b>	<b>4</b>
<b>4.1.</b>	<b><i>Theory of Operation.....</i></b>	<b>4</b>
<b>4.2.</b>	<b><i>Responsivity .....</i></b>	<b>4</b>
<b>4.3.</b>	<b><i>Dark Current.....</i></b>	<b>4</b>
<b>4.4.</b>	<b><i>Bandwidth and Response .....</i></b>	<b>5</b>
<b>4.5.</b>	<b><i>Terminating Resistance.....</i></b>	<b>5</b>
<b>4.6.</b>	<b><i>Gain Adjustment.....</i></b>	<b>5</b>
<b>Chapter 5</b>	<b>Troubleshooting .....</b>	<b>6</b>
<b>Chapter 6</b>	<b>Specifications .....</b>	<b>7</b>
<b>6.1.</b>	<b><i>Response Curve.....</i></b>	<b>9</b>
<b>6.2.</b>	<b><i>Mechanical Drawing .....</i></b>	<b>10</b>
<b>Chapter 7</b>	<b>Certificate of Conformance .....</b>	<b>11</b>
<b>Chapter 8</b>	<b>Regulatory.....</b>	<b>12</b>
<b>Chapter 9</b>	<b>Thorlabs Worldwide Contacts.....</b>	<b>13</b>

---

# Chapter 1 Warning Symbol Definitions

Below is a list of warning symbols you may encounter in this manual or on your device.

Symbol	Description
— — —	Direct Current
~	Alternating Current
~ ~	Both Direct and Alternating Current
⊕ ⊖	Earth Ground Terminal
○ ⊕	Protective Conductor Terminal
	Frame or Chassis Terminal
▽	Equipotentiality
	On (Supply)
○	Off (Supply)
■ ■	In Position of a Bi-Stable Push Control
■ ■ ■	Out Position of a Bi-Stable Push Control
	Caution: Risk of Electric Shock
	Caution: Hot Surface
	Caution: Risk of Danger
	Warning: Laser Radiation
	Caution: ESD Sensitive Components

## Chapter 2 Description

The PDA100A2 is an amplified, switchable-gain, Silicon (Si) detector designed for detection of light signals ranging from 320 nm to 1100 nm. An eight-position rotary switch allows the user to vary the gain in 10 dB steps. A buffered output drives 50 Ω load impedances up to 5 V. The PDA100A2 housing includes a removable threaded coupler (SM1T1) and retainer ring (SM1RR) that is compatible with any number of Thorlabs 1" threaded accessories. This allows convenient mounting of external optics, light filters, apertures, as well as providing an easy mounting mechanism using Thorlabs' cage assembly accessories. Also included is a ±12 V power supply.



### ESD Caution



The components inside this instrument are ESD sensitive. Take all appropriate precautions to discharge personnel and equipment before making any electrical connections to the unit.

## Chapter 3 Setup

The detector can be set up in many different ways using our extensive line of adapters. However, the detector should always be mounted and secured for best operation.

1. Unpack the optical head, install either an imperial or metric optical post into one of the mounting holes located on the bottom and side of the sensor, and mount using a post holder. Note that these detectors feature tapped holes that accept both 8-32 and M4 threads, so using either imperial or metric TR posts is possible.
2. Connect the power supply 3-pin plug into the power receptacle on the PDA100A2.
3. Plug the power supply into a 50 to 60Hz, 100 V / 120 V / 230 V power outlet.
4. Attach a  $50\ \Omega$  coax cable (i.e. RG-58U) to the output of the PDA. When running cable lengths longer than 12" we recommend terminating the opposite end of the coax with a  $50\ \Omega$  resistor (Thorlabs p/n T4119) for maximum performance. Connect the remaining end to a measurement device such as an oscilloscope or high speed DAQ card. **Caution:** Many high speed oscilloscopes have input impedances of  $50\ \Omega$ . In this case, do not install a  $50\ \Omega$  terminator. The combined loads will equal  $25\ \Omega$  which could allow  $\sim 135$  mA of output current. This will damage the output driver of the PDA100A2.
5. Power the PDA100A2 on using power switch located on top side of unit.

### Caution!


**The active detector area in the PDA100A2 is offset from the front of the housing. When using fiber adapters, make sure that the fiber ferrule does not crash into the detector. Failure to do so may cause damage to the diode and/or the fiber.**

6. Install any desired filters, optics, adapters, or fiber adapters to the input aperture.
7. During alignment, take appropriate precautions, such as using reduced radiation power, or other precautions, and use proper eye and/or skin protection as recommended by the radiation source manufacturer.
8. Apply a light source to the detector. Adjust the gain to the desired setting.

# Chapter 4 Operation

## 4.1. Theory of Operation

Thorlabs PDA series are ideal for measuring both pulsed and CW light sources. The PDA100A2 includes a reverse-biased PIN photo diode, mated to a switchable gain transimpedance amplifier, and packaged in a rugged housing.



## 4.2. Responsivity

The responsivity of a photodiode can be defined as a ratio of generated photocurrent ( $I_{PD}$ ) to the incident light power ( $P$ ) at a given wavelength:

$$R(\lambda) = \frac{I_{PD}}{P}$$

## 4.3. Dark Current

Dark current is leakage current which flows when a bias voltage is applied to a photodiode. The PDA with Transimpedance Amplifier does control the dark current flowing out. Looking at the figure above, it can be noted that Point B is held at ground and the amplifier will try to hold point A to "Virtual Ground". This minimizes the effects of dark current present in the system.

The dark current present is also affected by the photodiode material and the size of the active area. Silicon devices generally produce low dark current compared to germanium devices which have high dark currents. The table below lists several photodiode materials and their relative dark currents, speeds, sensitivity, and costs. Please note that sensitivity values in the table are typical values; Thorlabs offers photodetectors with sensitivity ranges that vary from those shown below.

Material	Dark Current	Speed	Sensitivity <sup>1</sup>	Cost
<b>Silicon (Si)</b>	Low	High	400 - 1000 nm	Low
<b>Germanium (Ge)</b>	High	Low	900 - 1600 nm	Low
<b>Gallium Phosphide (GaP)</b>	Low	High	150 - 550 nm	Med
<b>Indium Gallium Arsenide (InGaAs)</b>	Low	High	800 - 1800 nm	Med
<b>Extended Range: Indium Gallium Arsenide (InGaAs)</b>	High	High	1200 - 2600 nm	High

## 4.4. Bandwidth and Response

For the PDA detectors, the gain of the detector is dependent on the feedback element ( $R_f$ ). The bandwidth of the detector can be calculated using the following:

$$f(-3dB) = \sqrt{\frac{GBP}{4\pi R_f \times C_D}}$$

GBP is the amplifier gain bandwidth product and  $C_D$  is the sum of the photodiode junction capacitance and the amplifier capacitance.

## 4.5. Terminating Resistance

The maximum output of the PDA100A2 is 10 V for high impedance loads (i.e.  $R_{Load} > 5 \text{ k}\Omega$ ) and 5 V for 50  $\Omega$  loads. Adjust the gain so that the measured signal level out of the PDA100A2 is below 10 V (5 V with a 50  $\Omega$  load) to avoid saturation.

For low terminating resistors, <5 k $\Omega$  or 1% error, an additional factor needs to be considered. The output of the PDA includes a 50  $\Omega$  series resistor ( $R_s$ ). The output load creates a voltage divider with the 50  $\Omega$  series resistor as follows:

$$\text{Scale Factor} = \frac{R_{Load}}{R_{Load} + R_s}$$

$$V_{OUT} = \mathcal{R}(\lambda) * \text{Transimpedance Gain} * \text{Scale Factor} * \text{Input Power (W)}$$

Note that we already include the scale factor in our specification for the gain at 50  $\Omega$ . Refer to the table in Chapter 6 for additional performance specifications.

## 4.6. Gain Adjustment

The PDA100A2 includes a low noise, low offset, high gain transimpedance amplifier that allows gain adjustment over a 70 dB range. The gain is adjusted by rotating the gain control knob, located on the top side of the unit. There are 8 gain positions incremented in 10 dB steps. It is important to note that the bandwidth will decrease as the gain increases. See the specifications table on page 7 to choose the best gain versus bandwidth for a given input signal.

<sup>1</sup> Approximate values; actual wavelength values will vary.

## Chapter 5 Troubleshooting

Problem	Suggested Solutions
<b>There is no signal response.</b>	Verify that the power is switched on and all connections are secure.
	Verify the proper terminating resistor is installed if using a voltage measurement device.
	Verify that the optical signal wavelength is within the specified wavelength range.
	Verify that the optical signal is illuminating the detector active area.
<b>Output Voltage will not increase. Detector Output is skewed.</b>	Check to make sure the detector is not saturated. Refer to the Output Voltage spec. in the table on page 7.  Install a 1" Lens Tube (SM1L10) to the thread coupler (SM1T1) to baffle any external light sources to see if this improves the response.

# Chapter 6 Specifications

All performance specifications are typical, performed at 25 °C ambient temperature, and assume a 50 Ω load, unless stated otherwise.

Performance Specifications <sup>2</sup>			
0 dB Setting		40 dB Setting	
<b>Gain (Hi-Z)</b>	$1.51 \times 10^3$ V/A ±2%	<b>Gain (Hi-Z)</b>	$1.51 \times 10^5$ V/A ±2%
<b>Gain (50 Ω)</b>	$0.75 \times 10^3$ V/A ±2%	<b>Gain (50 Ω)</b>	$0.75 \times 10^5$ V/A ±2%
<b>Bandwidth<sup>3</sup></b>	11 MHz	<b>Bandwidth<sup>3</sup></b>	90 kHz
<b>Noise (RMS)</b>	268 μV	<b>Noise (RMS)</b>	229 μV
<b>NEP (@ λ<sub>p</sub>)</b>	$7.17 \times 10^{-11}$ W/√Hz	<b>NEP (@ λ<sub>p</sub>)</b>	$2.67 \times 10^{-12}$ W/√Hz
<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)	<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)
10 dB Setting		50 dB Setting	
<b>Gain (Hi-Z)</b>	$4.75 \times 10^3$ V/A ±2%	<b>Gain (Hi-Z)</b>	$4.75 \times 10^5$ V/A ±2%
<b>Gain (50 Ω)</b>	$2.38 \times 10^3$ V/A ±2%	<b>Gain (50 Ω)</b>	$2.38 \times 10^5$ V/A ±2%
<b>Bandwidth<sup>3</sup></b>	1.4 MHz	<b>Bandwidth<sup>3</sup></b>	28 kHz
<b>Noise (RMS)</b>	195 μV	<b>Noise (RMS)</b>	271 μV
<b>NEP (@ λ<sub>p</sub>)</b>	$6.75 \times 10^{-12}$ W/√Hz	<b>NEP (@ λ<sub>p</sub>)</b>	$4.2 \times 10^{-12}$ W/√Hz
<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)	<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)
20 dB Setting		60 dB Setting	
<b>Gain (Hi-Z)</b>	$1.5 \times 10^4$ V/A ±2%	<b>Gain (Hi-Z)</b>	$1.5 \times 10^6$ V/A ±5%
<b>Gain (50 Ω)</b>	$0.75 \times 10^4$ V/A ±2%	<b>Gain (50 Ω)</b>	$0.75 \times 10^6$ V/A ±5%
<b>Bandwidth<sup>3</sup></b>	800kHz	<b>Bandwidth<sup>3</sup></b>	9 kHz
<b>Noise (RMS)</b>	219 μV	<b>Noise (RMS)</b>	423 μV
<b>NEP (@ λ<sub>p</sub>)</b>	$3.36 \times 10^{-12}$ W/√Hz	<b>NEP (@ λ<sub>p</sub>)</b>	$6.24 \times 10^{-12}$ W/√Hz
<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)	<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)
30 dB Setting		70 dB Setting	
<b>Gain (Hi-Z)</b>	$4.75 \times 10^4$ V/A ±2%	<b>Gain (Hi-Z)</b>	$4.75 \times 10^6$ V/A ±5%
<b>Gain (50 Ω)</b>	$2.38 \times 10^4$ V/A ±2%	<b>Gain (50 Ω)</b>	$2.38 \times 10^6$ V/A ±5%
<b>Bandwidth<sup>3</sup></b>	260 kHz	<b>Bandwidth<sup>3</sup></b>	3 kHz
<b>Noise (RMS)</b>	222 μV	<b>Noise (RMS)</b>	1.22 mV
<b>NEP (@ λ<sub>p</sub>)</b>	$2.83 \times 10^{-12}$ W/√Hz	<b>NEP (@ λ<sub>p</sub>)</b>	$7.88 \times 10^{-12}$ W/√Hz
<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)	<b>Offset</b>	±8 mV (Typ.) ±12 mV (Max)

<sup>2</sup> The PDA100A2 has a 50 Ω series terminator resistor (i.e. in series with amplifier output). This forms a voltage divider with any load impedance (e.g. 50 Ω load divides signal in half).

<sup>3</sup> Tested at 650 nm wavelength. For NIR wavelengths, the rise time of the photodiode element will become slower which may limit the effective bandwidth of the amplified detector.


Electrical Specifications		
<b>Detector</b>	-	Si
<b>Active Area</b>	-	Ø9.8 mm (75.4 mm <sup>2</sup> )
<b>Wavelength Range</b>	$\lambda$	320 to 1100 nm
<b>Peak Wavelength</b>	$\lambda_p$	960 nm (Typ.)
<b>Peak Response</b>	$R(\lambda_p)$	0.72 A/W (Typ.)
<b>Amplifier GBP</b>	-	600 MHz
<b>Output Impedance</b>	-	50 Ω
<b>Max Output Current</b>	I <sub>OUT</sub>	100 mA
<b>Load Impedance</b>	-	50 Ω to Hi-Z
<b>Gain Adjustment Range</b>	-	0 to 70 dB
<b>Gain Steps</b>	-	8 x 10 dB Steps
<b>Output Voltage</b>	V <sub>OUT</sub>	0 to 5 V (50 Ω) 0 to 10 V (Hi-Z)

General	
<b>On/Off Switch</b>	Slide
<b>Gain Switch</b>	8 Position Rotary
<b>Output</b>	BNC (DC Coupled)
<b>Package Size</b>	2.79" x 2.07" x 0.89" (70.9 mm x 52.5 mm x 22.5 mm)
<b>PD Surface Depth<sup>4</sup></b>	0.13" (3.3 mm)
<b>Weight, Detector Only</b>	0.10 kg
<b>Accessories</b>	SM1T1 Coupler SM1RR Retainer Ring
<b>Operating Temp</b>	10 to 40 °C
<b>Storage Temp</b>	-20 to 70 °C
<b>AC Power Supply</b>	AC - DC Converter
<b>Input Power<sup>5</sup></b>	6 W 100 V / 120 V / 230 V, 50 - 60 Hz


<sup>4</sup> Measured from the front of the housing body.

<sup>5</sup> Although the power supply is rated for 6 W, the detector's actual usage is <5 W over the full operating range.

## 6.1. Response Curve



## 6.2. Mechanical Drawing



## Chapter 7 Certificate of Conformance



# THORLABS

[www.thorlabs.com](http://www.thorlabs.com)

### EU Declaration of Conformity

*in accordance with EN ISO 17050-1:2010*

We: Thorlabs Inc.

Of: 56 Sparta Avenue, Newton, New Jersey, 07860, USA

*in accordance with the following Directive(s):*

2014/35/EU Low Voltage Directive (LVD)

2014/30/EU Electromagnetic Compatibility (EMC) Directive

2011/65/EU Restriction of Use of Certain Hazardous Substances (RoHS)

*hereby declare that:*

**Model:** PDA05CF2, PDA100A2, PDA10A2, PDA10CS2, PDA10D2, PDA20CS2, PDA25K2, PDA30B2, PDA36A2, PDA50B2, PDA8A2, PDA20C2

**Equipment:** *Switchable adjustable and Non-switchable fixed gain Amplified Photodetector*

*is in conformity with the applicable requirements of the following documents:*

EN 61010-1	Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use.	2010
------------	---	------

EN 61326-1	Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements	2013
------------	---	------

*and which, issued under the sole responsibility of Thorlabs, is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8th June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, for the reason stated below:*

*does not contain substances in excess of the maximum concentration values tolerated by weight in homogenous materials as listed in Annex II of the Directive*

*I hereby declare that the equipment named has been designed to comply with the relevant sections of the above referenced specifications, and complies with all applicable Essential Requirements of the Directives.*

*Signed:*

*On:* 09 January 2018

Name: Ann Strachan

Position: Compliance Manager

EDC - PDA05CF2, PDA100A2, PDA10A2, P...




## Chapter 8 Regulatory

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return "end of life" units without incurring disposal charges.

- This offer is valid for Thorlabs electrical and electronic equipment:
- Sold after August 13, 2005
- Marked correspondingly with the crossed out "wheelie bin" logo (see right)
- Sold to a company or institute within the EC
- Currently owned by a company or institute within the EC
- Still complete, not disassembled and not contaminated

As the WEEE directive applies to self contained operational electrical and electronic products, this end of life take back service does not refer to other Thorlabs products, such as:



**Wheelie Bin Logo**

- Pure OEM products, that means assemblies to be built into a unit by the user (e.g. OEM laser driver cards)
- Components
- Mechanics and optics
- Left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.

### **Waste Treatment is Your Own Responsibility**

If you do not return an "end of life" unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

### **Ecological Background**

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of life products will thereby avoid negative impacts on the environment.

# Chapter 9    Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at  
[www.thorlabs.com/contact](http://www.thorlabs.com/contact) for our most up-to-date contact information.



## **USA, Canada, and South America**

Thorlabs, Inc.  
[sales@thorlabs.com](mailto:sales@thorlabs.com)  
[techsupport@thorlabs.com](mailto:techsupport@thorlabs.com)

## **UK and Ireland**

Thorlabs Ltd.  
[sales.uk@thorlabs.com](mailto:sales.uk@thorlabs.com)  
[techsupport.uk@thorlabs.com](mailto:techsupport.uk@thorlabs.com)

## **Europe**

Thorlabs GmbH  
[europe@thorlabs.com](mailto:europe@thorlabs.com)

## **Scandinavia**

Thorlabs Sweden AB  
[scandinavia@thorlabs.com](mailto:scandinavia@thorlabs.com)

## **France**

Thorlabs SAS  
[sales.fr@thorlabs.com](mailto:sales.fr@thorlabs.com)

## **Brazil**


Thorlabs Vendas de Fotônicos Ltda.  
[brasil@thorlabs.com](mailto:brasil@thorlabs.com)

## **Japan**

Thorlabs Japan, Inc.  
[sales@thorlabs.jp](mailto:sales@thorlabs.jp)

## **China**

Thorlabs China  
[chinasales@thorlabs.com](mailto:chinasales@thorlabs.com)



**THORLABS**  
[www.thorlabs.com](http://www.thorlabs.com)

