


## Product Presentation



## Application

GT-N19 series Paper and cardboard thickness is measured the thickness of paper and cardboard special apparatus, widely used in paper, cardboard and other sheet materials, measured the thickness of corrugated cardboard.

The instrument is the cardboard is the cardboard carton production, research and inspection of enterprises and sectors, such as an indispensable instrument used.

## Standards

ISO 534, ISO 3034, GB/T 451.3, GB/T 6547

## Key Specification

Item	GT-N19A	GT-N19B	GT-N19C	GT-N19D
Measuring range	0-6mm			
Accuracy	0.01mm		0.001mm	
Contact pressure	100±1kPa			
Contact area	200±5mm2			
Dial	Pointer Dial indicator	Digital Dial indicator	Pointer Dial gauge	Digital Dial gauge
Parallelism	≤0.002mm			
Dimension	375x264x400			
Weight	About 6kg			