

## 90-1 Centrifuge

Outline of uses:

90-1 is an improvement on 800. It is possessed of these special features following:

It is beautiful in appearance. Large in capacity but small in bulk and has a complete function and a stable performance, you can chose the speed that you need at the set range, and it can keep balance automatically, in addition, it has a little rate of increase in temperature and can be used efficiently and widely.

It is a dependable instrument for qualitative analysis of serum, plasma and immunity in hospitals or chemical and biochemical laboratories.



| Model | MaxSpeed  | Capacity | MaxRCF | TimerRange | OutlineDimension | NetWeight |
|-------|-----------|----------|--------|------------|------------------|-----------|
| 90-1  | 4000r/min | 20ml*6   | 1795g  | 0-60min    | 280*240*260mm    | 7Kg       |