

EN 71 Part 1 Length of Ropes and Chains Tester GT-MB22A

Product Introduction

Application:

Used to measure length for the ropes, chains or wires on toys

Standards:

EN-71-1 section 8.20 , 8.40

ISO 8124-1 section 5.11.1

GB6675-2 section 5.11.1

Dimensions

External Dimension:	
500*500*175 mm(L*W*H)	
Package Dimension:	
Gross Weight:	
Net Weight:	
16KG	

Key Specifications

Model	GT-MB22	
Scale	1500mm	
Test load	25N	
Dimensions(L x D x H)	50x50x175cm	
Weight	16kg	
Standards	EN-71-1 section 8.20 , 8.40 ISO 8124-1 section 5.11.1 GB6675-2 section 5.11.1	
Standard accessory	25N weight	1pc