

Colonometer

>>> J-2 & J-3 Colonometer

Features:

- Compact structure makes it less space occupation. Adapted for various environment.
- The LED light can be adjusted to make an ideal perspective for the operator.
- It will beeps while counting, which can avoid wrong counting.
- Unique Up-light source makes it convenient to count bacteria that is hard to observe.
- Specified designed counting pen to ensure no omit of the cells.
- Wide-range voltage design. Avoid the interference to the equipment by instable voltage.

Technical Data:

	J-2	J-3
Lighting:	Matrix LED(White light)	Matrix LED(White light)
Digital display:	3 digits	3 digits
Counting range:	0-999	0-999
Compatible petri dish:	50-90mm	50-150mm
Input voltage(frequency):	AC 100-240V (50/60Hz)	AC 100-240V (50/60Hz)
Input power:	20W	40W
Magnification:	3-9 Times	3-9 Times
Protective class:	IP21	IP21
Permissible relative moisture:	80%	80%
Permissible ambient temperature:	5-50°C	5-50°C
Dimensions:	255×210×160mm	360×300×180mm
Net. Weight:	2.2Kg	4.0Kg

Product-details:

