

WX Horizontal Autoclave



Features:

- ◆ Full SUS 304 stainless steel structure
- ◆ The control system adopts imported original, Chamber pressure and inter-layer pressure can be monitored
- ◆ The vacuum pump is imported from Siemens, make the sterilization effect more thoroughly.
- ◆ Automatic quickly open the door, safety interlock system.
- ◆ LCD dynamic display temperature and time
- ◆ Error analysis, fault warning display system
- ◆ Two pressure gauge indicating inner and interlayer pressure
- ◆ With a variety of optional sterilization procedures
- ◆ Automatic overpressure release and power off protection, safety ensuring at all time

| Model | WXGJDS-0.4A | WXGJDS-0.6A | WXGJDS-0.8A | WXGJDS-1.0A | WXGJDS-1.2A |
|-------------------------|-------------|-------------|-------------|-------------|-------------|
| Volume | 400L | 600L | 800L | 1000L | 1200L |
| Max Working Temperature | 134°C | 134°C | 134°C | 134°C | 134°C |

| | | | | | |
|-----------------------|----------------|-----------------|-----------------|-----------------|------------------|
| Max Working Pressure | 0.22Mpa | 0.22Mpa | 0.22Mpa | 0.22Mpa | 0.22Mpa |
| Sterilization Chamber | 962*600*700 mm | 1430*600*700 mm | 1850*600*700 mm | 1700*680*910 mm | 1700*680*1050 mm |
| G.W/N.W | 1050kg/950kg | 1100kg/1000kg | 1300kg/1200kg | 1500kg/1400kg | 2080kg/1970kg |
| Power | 380v 20kW | 380v 45kW | 380v 45kW | 380v 50kW | 380v 50kW |