



U-MEDICAL
合众物联

BSC-18L/28L

Portable/Vehicle refrigerator





U-MEDICAL
合众物联

Portable Vehicle refrigerator

Portable lightweight design

- 1. ABS Engineering Plastic Shell, 304 stainless steel lining, optional built-in wire basket, easy classification storage.
- 2. Standard Aluminum Alloy handle, strap fastening, shoulder back, hand can be easy to be fixed in the car trunk. (BSC-18L)

German technology, Quality of craftsmanship

- 1. AC POWER SUPPLY VOLTAGE AC 220V ~ , Vehicle power supply voltage DC12/24V, perfect battery protection.
- 2. Germany SECOP DC compressor, Germany EBM Fan, micro channel condenser, high heat transfer efficiency, stable operation, especially suitable for vehicle transport environment.
- 3. The air-cooled circulating fan is arranged in the box, the temperature is stable and uniform, and the fluctuation is small, Built-in LED lighting , easy to operate.
- 4. High-pressure foaming insulation layer, Non-fluorine refrigerant.
- 5. Temperature setting range:-20 °C ~ 10 °C, accuracy 0.1 °C (ambient temperature 32 °c).

High precision control system, more security guarantee

- 1. Embedded intelligent control system, 2.8-inch TFT touch screen, easy to operate, dynamic display running status.
- 2. Optional 4G wireless networking function, Docking cold-chain monitoring cloud platform, storage temperature monitoring anytime, anywhere, one more security protection.
- 3. Optional Bluetooth Printer, easy to print temperature history at any time.
- 4. Support USB Data Export Function, can export PDF, excel tables, support the generation of statistical reports or charts.

Model	External size (mm)	Capacity (L)	Temperature (°C)	Power (W)	Net Weight (kg)
BSC-18L	500*350*350	18	10 ~ -20	42	10.5
BSC-28L	500*350*440	28	10 ~ -20	50	12